



# 台灣日研股份有限公司

338桃園市蘆竹區中興一街67號

## 產品規格表

品名:鰹魚抽出粉T-28

項目/Item	規格/Specification	分析方式/Analysis <sup>1</sup> Method	
水分/Moisture	6.0% max	Loss on Drying	
鹽分/Salt	15.0%±3.0%	In house method modified from AOAC	
總氮數/Total Plate Count	2.4%±0.5%	Kjeldahl method (based on AOAC)	
總落菌數/Total Plate Count	3,000 cfu/g max	Based on FDA BAM online	
大腸桿菌/Coliforms	Negative	Based on FDA BAM online	
黴菌&酵母菌Yeast&Mold	200 cfu/g max	Based on FDA BAM online	
沙門氏菌/Salmonella	Not detected	Based on FDA BAM online	
金黃色葡萄球菌S.aureus	100 cfu/g max	Based on FDA BAM online	
霍亂弧菌V.cholerae	Not detected	Based on FDA BAM online	
鎘Cadmium(ppm)	≤0.1	In house method based on AOAC	
鉛/Lead(ppm)	≤0.3	In house method based on AOAC	
無機砷Inorganic Arsenic(ppm)	≤2.0	In house method based on AOAC	
汞Mercury(ppm)	≤0.4	In house method based on Food USEPA and ASTM	
錫Tin(ppm)	≤250	In house method based on AOAC	
鋅Zinc(ppm)	≤100	In house method based on AOAC	
銅Copper(ppm)	≤20.0	In house method based on AOAC	
黃麴毒素Aflatoxin(B )(ppb)	≤10	In house method SOP based on QuEChERS method	
苯甲酸(Benzoic acid)	≤200	In house method SOP based on QuEChERS method	
組織胺(Histamine)	≤2.0	In house method SOP based on QuEChERS method	
水活性(Water Activity)	≤0.7	In house method based on AOAC	
感官測試/Sensory Test	Proper to key sample	Sensory Test	
淨重: 18公斤 (18公斤x1)		保存期限: 在未開封並妥善保存之情況下為18個月。	
貯存條件：需密封，保存於陰涼，避免光照並乾燥的地方。			
生產程序：			
(生料)→ 檢查 → 混合加熱(≥92℃,≥2min)→ 調整 →過篩（150目）→ 磁鐵 <sup>1</sup> → 乾燥噴霧 → 過篩（50目）			
磁鐵 <sup>1</sup> → 檢查 → 包裝 → 檢查 → 出貨			
	組成原料	原料原產地	過敏原 <sup>3</sup>
麥芽糊精 (Maltodextrin)	木薯	泰國	
鰹魚抽出物 (Bonito Extract)	鰹魚	泰國	None
鹽 Salt	岩鹽	泰國	
柴魚乾(Dried Bonito Extract)	鰹魚	泰國	
酵母抽出	酵母抽出 岩鹽	泰國 泰國	
備註:/Remark：			
**：本版規格書係依泰國日研原廠英文版規格翻譯,中文譯本僅供參考,文義如與英文有歧異,概以英文版規格書為準。			

本規格資訊以發布日期為準，並可能隨時進行更改修正。

\* 1：可以使用對等性檢驗方法。

\* 2：本產品每年會經過研發部分析一次。成分依其含量多寡由高至低標示。

\* 3：成分依其含量多寡由高至低標示。

\* 4：過敏原聲明是依據台灣食品過敏原標示規定（依據食品安全衛生管理法第22條第1款第10項）。

※ 括號中之過敏原為非強制性。